

CMPRK


CM ProLine



Article nr: CMPRK

- ✓ HYGIENIC
- ✓ COMFORTABLE
- ✓ STABLE

The backcushion can be used to prevent the patient from descending in the mattress when x-ray is being used from the side. Two handles facilitates handling the cushion.



PRODUCTINFORMATION:

Art.nr: CMPRK
Coating: Artificial leather
(Technoskin)
Core: Polyether/EVA
Two handles.
Size: 55x36 cm Height: 16-8 cm

Cleaning: Disinfection with
alcohol/soap solution.



Distributor:

COMFORT medical

Manufactured by:

Comfort Medical AB
www.comfortmedical.se