


# HIGH POSITIONING-CUSHIONS

CM ProLine



The postoperative highpositioning-cushion has a bowl-shaped top. Can also be made with a plane top.



The peroperative highpositioning-cushion has a core of sensaten wich makes the pillow hard and stable.

- ✓ HYGIENIC
- ✓ COMFORTABLE
- ✓ STABLE

## PRODUCTINFORMATION:

Post operative cushion

Coating: Technoskin (artificial leather).

Core: Polyether.

1. Size: 80x30x22 cm bowl-shaped

Article nr: CMPORTO80

2. Size: 60x30x22 cm

Article nr: CMPORTO60

Pre operative cushion

Coating: Technoskin (artificial leather).

Core: Sensaten (hard).

1. Size: 80x22x22 cm (plane) .

Article no: CMPORTO80H

2. Size: 60x22x22 cm

Article no: CMPORTO60H

**Cleaning:** Disinfection with alcohol/ soap solution.



**DISTRUBUTOR:**

**COMFORT** medical

**Manufactured by:**

Comfort Medical AB

[www.comfortmedical.se](http://www.comfortmedical.se)