

Knee cushions

CM ProLine



Article nr: CMPKK
Name: Kneecushion
Size: 45x31cm Height: 17 cm

- ✓ Hygienic
- ✓ Comfortable
- ✓ Ergonomic

PRODUCT INFORMATION:

Kneecushions:

Coating: Technoskin artificial leather

Core: Polyether

Two handles on cushion with technoskin.
Can be made with velcro sealed solution for MR-environment.


Small Kneecushions:

Has welded PU-coating. Core of polyether (soft) or sensaten (hard). The 40 cm long version can also be used for CT-colon in prone position to avoid compression of colon transversum.

Cleaning: Disinfection with alcohol/soap solution.



Article nr: CMPKKF
Name: Kneecushion Flat
Size: 55x38 cm
Height: 16-8 cm



Article nr: CMPKKC
Name: Cylindrical Kneecushion
Size: 50x15 cm



Article nr: CMPLKM25 (soft)
Name: Small Kneecushion 25
Size: 25x17 cm
Height: 8 cm

Article nr: CMPLKM40 (soft)
CMPLKH40 (hard)
Name: Small Kneecushion 40
Size 40x17 cm
Height 8 cm

Distributor:

Manufactured by:
Comfort Medical AB
www.comfortmedical.se

COMFORT medical