

POSITIONING CUSHIONS

CM ProLine



- ✓ Hygienic
- ✓ Optimum grip
- ✓ Gives no artifacts

PRODUCT INFORMATION:

The cushions are made of homogenous EVA and sensaten with closed cells.


High friction material, gives no artifact. CM ProLine comes with 1 year function guarantee.

Care:

The material does not absorb any fluid and is easy to rinse. CM ProLine is also easy to clean with soap-alcohol solutions.

Article no:


See on page 2.




Comfort Medical AB
Hjalmar Brantings väg 11G
SE - 261 43 Landskrona
Sweden

Mail: info@comfortmedical.se
www.comfortmedical.se


COMFORT medical




Art.nr: CMP60 Wedge 60 x 25 x 12 cm, 45°




Art.nr: CMP50 Wedge 50 x 25 x 12 cm, 45°




Art.nr: CMP40 Wedge 40 x 25 x 12 cm, 45°




Art.nr: CMPF40 Wedge flat 40 x 22 x 12 cm, 30°




Art.nr: CMP2040 block with slot for cassetes. 40 x 20 x 5 cm




Art.nr: CMPHK Handcushion 18 x 15 x 6 cm



Art.nr: CMPHT Stair shaped
fingercushion 13 x 13 x 8 cm



Art.nr: CMPKS wedge low
19 x 19 x 0,5-4 cm



Art.nr: CMPKL wedge high
19 x 19 x 0,5-6 cm



Art.nr: CMPHLK Handcushion
23 x 16 x 10 cm, 45° angle