

CM X-Ray Mattress

CM ProLine



- ✓ Comfortable
- ✓ Low attenuation
- ✓ Durable
- ✓ Ergonomic

The X-ray Mattres can be manufactured in different shapes and sizes. Low friction underneath and eight handles for easy adjustment of the mattres.



PRODUCT INFORMATION:

Coating: Technoskin - artificial leather.

Core: Viscoelastic polyurethanfoam/EVA sandwich.

Eight handles, low friction underneath.

Thickness 3 cm.

Cleaning: Disinfection with alcohol/soap solution.

Article nr: CM + width and length
For example: CM65220



We deliver mattresses to several x-ray companies.

Distributor:

COMFORT medical

Manufactured by:
Comfort Medical AB
www.comfortmedical.se