

MAVIG - mobile, height adjustable X-Ray protective bed- shield

- ▶ movable, tilt resistance design, absolutely smooth running and flexible
- ▶ special X-base with plastic covers for minimum required space
- ▶ X-Ray protective material on natural rubber base with double-sided, kind to the skin and low-maintenance surface coating
- ▶ ideal adjustment for the respective requirements on site through the user-friendly adjustability of the height of the shield (sliding lift through the pressure spring assistance)



In accordance with the needs for minimization ever more frequently there are requirements for practicable radiation protection screens at the patient-beds.

A goal-oriented solution, to be integrated easily into the usual work routine and which considers quite possible spatially close local dimensionings, is presented with this new development of our company.

Already with a L.E. value of normally 0.5 mm an effective radiating reduction is obtained. The dimensions of the screen can be co-ordinated individually with the requirements of the specific application.



MAVIG GmbH
Strahlenschutz und
Systemtechnik

Postal address: P.
O. Box 82 03 62
81803 Munich
Germany
Address:
Stahlgruberring 5
81829 Munich
Germany

Phone:
0 89 / 4 20 96 -0

Fax:
0 89 / 4 20 96 -200

eMail:
info@mavig.com
Commercial
register
Munich
HRB-Nr. 76 403

tax-no.:
817/27263

Manager:
Manfred Stoian
Christian Stoian

VAT-UD No:
DE811268845